

Oxford Accupet EVO Single Channel Pipette

Volume lock setting to prevent accidental volume drifts

Made using highly durable material for corrosion and physical shock resistance

Highly durable PVDF tip cones



Corrosion resistant tip ejector with a unique shock absorbing mechanism



Magnetically assisted piston for precise and consistent results



Innovative spring and seal design for one of the lowest plunger forces for a fatigue free experience



Universal tip compatability

Model#	Description	Volume	No. of Channels
AE-2.5	Accupet EVO Single Channel Adjustable Vol. Pipette	0.1-2.5 μl	Single
AE-10	Accupet EVO Single Channel Adjustable Vol. Pipette	0.5-10 μΙ	Single
AE-20	Accupet EVO Single Channel Adjustable Vol. Pipette	2-20 µl	Single
AE-50	Accupet EVO Single Channel Adjustable Vol. Pipette	5-50 μΙ	Single
AE-100	Accupet EVO Single Channel Adjustable Vol. Pipette	10-100 μΙ	Single
AE-200	Accupet EVO Single Channel Adjustable Vol. Pipette	20-200 μΙ	Single
AE-1000	Accupet EVO Single Channel Adjustable Vol. Pipette	100-1000 μΙ	Single
AE-5000	Accupet EVO Single Channel Adjustable Vol. Pipette	500-5000 μΙ	Single
AE-10000	Accupet EVO Single Channel Adjustable Vol. Pipette	1000-10000 μΙ	Single

3 Year Warranty

For more details, visit us at www.oxfordlp.com



Oxford Accupet EVO Multichannel Pipette



30 - 300µl



Spring-loaded tip cone for less strain & proper tip fitment for precise pipetting



Color coding for easy identification



Large and clear 4-digit display

Hiahlv	durable	PVDF tir	cones

Model#	Description	Volume	No. of Channels
AE8-10	Accupet EVO 8-Channel Adjustable Vol. Pipette	0.5-10 μΙ	8
AE8-20	Accupet EVO 8-Channel Adjustable Vol. Pipette	2-20 μΙ	8
AE8-50	Accupet EVO 8-Channel Adjustable Vol. Pipette	5-50 μl	8
AE8-100	Accupet EVO 8-Channel Adjustable Vol. Pipette	10-100 μΙ	8
AE8-200	Accupet EVO 8-Channel Adjustable Vol. Pipette	20-200 μΙ	8
AE8-300	Accupet EVO 8-Channel Adjustable Vol. Pipette	30-200 μΙ	8
AE12-10	Accupet EVO 12-Channel Adjustable Vol. Pipette	0.5-10 μΙ	12
AE12-20	Accupet EVO 12-Channel Adjustable Vol. Pipette	2-20 μΙ	12
AE12-50	Accupet EVO 12-Channel Adjustable Vol. Pipette	5-50 μΙ	12
AE12-100	Accupet EVO 12-Channel Adjustable Vol. Pipette	10-100 μΙ	12
AE12-200	Accupet EVO 12-Channel Adjustable Vol. Pipette	20-200 μΙ	12
AE12-300	Accupet EVO 12-Channel Adjustable Vol. Pipette	30-300 μl	12

3 Year Warranty

For more details, visit us at www.oxfordlp.com