


Accupet Pro

Single Channel Pipette

Designed to offer optimal performance and precision with effortless handling. The Accupet Pro series of pipettes are ideal for long hours of pipetting. It's intelligent design includes features like one-handed tip ejection and lightweight construction, offering exceptional comfort and convenience.



AccuPet Pro Single Channel Pipette Specifications

Single Channel Catalog No.	Volume Range (µL)	Increments (µL)	Test Volume (µL)	Inaccuracy (±) %	Imprecision (±) %
AP-2.5	0.1 - 2.5	0.01	0.25 2.5	12 2.5	6 1.6
AP-10	0.5-10	0.1	1 10	2.5 1	1.5 .08
AP-20	2 - 20	0.1	2 20	3 .09	2 0.4
AP-50	5 -50	0.5	5 50	2 0.6	2 0.3
AP-100	10 to 100	0.5	10 100	3 0.8	1.5 0.15
AP-200	20-200	1	20 200	2 0.6	0.8 0.15
AP-1000	100-1000	5	100 1000	2 0.6	0.7 0.2
AP-5000	500-5000	50	500 5000	2 0.5	0.6 0.2
AP-10ML	500-10000	100	500 10000	3 0.6	0.6 0.2

Also Available In Multichannel

For more details, visit us at www.oxfordlp.com