


Accupet Pro Multichannel Pipette

Designed to offer optimal performance and precision with effortless handling. The Accupet Pro series of pipettes are ideal for long hours of pipetting. Its intelligent design includes features like one-handed tip ejection and lightweight construction, offering exceptional comfort and convenience.



AccuPet Pro Multichannel Pipette Specifications

8-Channel Catalog No.	12-Channel Catalog No.	Volume Range (µL)	Increments (µL)	Test Volume (µL)	Inaccuracy (±) %	Imprecision (±) %
AP8-10	AP12-10	0.5 to 10	0.1	1 10	4 1.5	4 1.5
AP8-50	AP12-50	5 to 50	0.5	5 50	3 1	2 0.7
AP8-100	AP12-100	10 to 100	1	10 100	1.5 1	1.5 0.5
AP8-200	AP12-200	20 to 200	1	20 200	1.5 0.7	0.75 0.25
AP8-300	AP12-300	30 to 300	1	30 300	1.5 0.8	0.75 0.25

Also Available In Single Channel

For more details, visit us at www.oxfordlp.com