

Oxford BenchMate Plus

Single-Channel Ultra Smooth Pipettes

Classic Quality with a Modern Design

The BenchMate pipette line has been synonymous with durability and quality for over 69 years. The Oxford BenchMate Plus Pipette is a modern interpretation of the durable design, adjusting the material and body to create a lighter and more comfortable pipetting experience. Despite these updates, the BenchMate Plus retains the critical durability and inherent quality that users recognize as soon as they pick it up.

Incredible Smoothness and Precision

Any pipette can be precise, however a truly special pipette can help the user be more precise. The BenchMate Plus offers an incredibly smooth piston design assisted by an internal magnetic system. This design allows the plunger to be operated with minimal resistance and the magnetic system helps users to identify the first stop and second stop, ensuring that every pipetting stroke is precise beyond just the mechanical specifications. Users can pipette confidently and comfortably at every use.



SKU	Description	Volume Range
OBP-2.5	Oxford BenchMate Plus Single Channel Pipette. 0.1-2.5 µl	0.1-2.5 µl
OBP-10	Oxford BenchMate Plus Single Channel Pipette. 0.5-10 µl	0.5-10 µl
OBP-20	Oxford BenchMate Plus Single Channel Pipette. 2-20 µl	2-20 µl
OBP-50	Oxford BenchMate Plus Single Channel Pipette. 5-50 µl	5-50 µl
OBP-100	Oxford BenchMate Plus Single Channel Pipette. 10-100 µl	10-100 µl
OBP-200	Oxford BenchMate Plus Single Channel Pipette. 20-200 µl	20-200 µl
OBP-1000	Oxford BenchMate Plus Single Channel Pipette. 100-1000 µl	100-1000 µl
OBP-5000	Oxford BenchMate Plus Single Channel Pipette. 500-5000 µl	500-5000 µl
OBP-10000	Oxford BenchMate Plus Single Channel Pipette. 1000-10000 µl	1000-10000 µl

Oxford BenchMate Plus Multichannel Pipettes

Multi-Channel Ultra Smooth Pipettes

Precision and Comfort

Oxford BenchMate Plus Multichannel Pipettes offer an incredibly smooth pipetting experience typically missing in the multichannel pipette design. The mechanical specifications inherently offer advanced accuracy and precision, but the smooth piston design adds an additional layer of precision. The BenchMate Plus Multichannel Pipette is designed so that users can easily pipette 8 or 12 channels without excessive tip attachment or plunging force. The soft plunging force means that users can repeatedly and consistently pipette without combatting the typical multichannel resistance.

Efficiency Made Easy

Multichannel pipettes are good for reducing the workload needed to do a large amount of pipetting, however the typical multichannel pipette requires more force and is less comfortable to use. The Oxford BenchMate Plus Multichannel Pipette tackles this exact issue by reducing the plunger force and shaping the pipette to fit the user's hand. This premium pipette offers comfort and efficiency without sacrificing accuracy or precision.



SKU	Description	Volume Range
OBP8-10	Oxford BenchMate Plus Multi-Channel Pipette. 8-channel, 0.5-10 µl	0.5-10 µl
OBP8-20	Oxford BenchMate Plus Multi-Channel Pipette. 8-channel, 2-20 µl	2-20 µl
OBP8-50	Oxford BenchMate Plus Multi-Channel Pipette. 8-channel, 5-50 µl	5-50 µl
OBP8-100	Oxford BenchMate Plus Multi-Channel Pipette. 8-channel, 10-100 µl	10-100 µl
OBP8-200	Oxford BenchMate Plus Multi-Channel Pipette. 8-channel, 20-200 µl	20-200 µl
OBP8-300	Oxford BenchMate Plus Multi-Channel Pipette. 8-channel, 30-300 µl	30-300 µl

SKU	Description	Volume Range
OBP12-10	Oxford BenchMate Plus Multi-Channel Pipette. 12-channel, 0.5-10 µl	0.5-10 µl
OBP12-20	Oxford BenchMate Plus Multi-Channel Pipette. 12-channel, 2-20 µl	2-20 µl
OBP12-50	Oxford BenchMate Plus Multi-Channel Pipette. 12-channel, 5-50 µl	5-50 µl
OBP12-100	Oxford BenchMate Plus Multi-Channel Pipette. 12-channel, 10-100 µl	10-100 µl
OBP12-200	Oxford BenchMate Plus Multi-Channel Pipette. 12-channel, 20-200 µl	20-200 µl
OBP12-300	Oxford BenchMate Plus Multi-Channel Pipette. 12-channel, 30-300 µl	30-300 µl