

Oxford LTR Pipette Tips Reload Stacks

LTS[™]-Compatible Pipette Tips Refill Wafers

Premium Pipetting, Sustainable Packaging

Oxford LTR Pipette Tips Reload Stacks offer LTS[™] level performance in a sustainable packaging and more accessible price. Save money and the environment with refillable tip wafers stacked into a tower. Existing pipette tip racks can be utilized with refill wafers by removing the original wafer and replacing it with a refill wafer.

Optimal Sample Delivery

Every Oxford LTR Pipette Tip features a hydrophobic interior and an ultra-fine tip point that maximizes sample delivery. Furthermore, all Oxford LTR Pipette Tips are certified free of RNase, DNase, Cytotoxins, PCR inhibitors, endotoxins, and pyrogens.

High Quality Manufacturing

Our LTR Pipette Tips are manufactured in a fully automatic, touch-free facility using diamond-finished molds and high-quality medical grade polypropylene. Various QC tests and regulations ensure that each tip is manufactured according to rigorous Oxford Lab Products high quality standards.



Volumes				
Cat. No.	Description	Tips/Reload	Reloads/ Stack	Stacks/ Case
LTR-20-LR	20uL LTS-Compatible, Low-Retention Pipette Tips Reload Stack.96 Tips/Reload, 10 Reloads/Stack, 10 Stacks/Case	96	10	10
LTR-200-LR	200uL LTS-Compatible, Low-Retention Pipette Tips Reload Stack. 96 Tips/Reload, 10 Reloads/Stack, 10 Stacks/Case	96	10	10
LTR-1000-LR	1000uL LTS Compatible, Low Retention Pipette Tips Reload Stack. 96 Tips/Reload, 5 Reloads/Stack, 10 Stacks/Case	96	5	10



Product appearance, catalog numbers, prices, specifications, and technical information are subject to change without notice. oxfordlp.com • info@oxfordlp.com • 1-800-554-4920