

Oxford LTR Pipette Tips

LTS™-Compatible Pipette Tips

Performance You Expect, Price You Want

Oxford LTR Pipette Tips offer LTS™ level performance at a much more affordable price. These tips are engineered to securely attach and seal onto LTS™/XLS pipettes. Every tip features an ultra-fine tip point that guarantees accuracy, precision, and consistency.

Optimal Sample Delivery

Every tip features a hydrophobic interior to ensure that minimal sample is retained, and the full volume is delivered accurately. Furthermore, all Oxford LTR Pipette Tips are certified free of RNase, DNase, Cytotoxins, PCR inhibitors, endotoxins, and pyrogens.

High Quality Manufacturing

Our LTR Pipette Tips are manufactured in a fully automatic, touch-free facility using diamond-finished molds and high-quality medical grade polypropylene. Various QC tests and regulations ensure that each tip is manufactured according to rigorous Oxford Lab Products high quality standards.



Volumes

Cat. No.	Description	Tips/Rack	Racks/Pack	Packs/Case
LTR-20-L	20uL LTS-Compatible Low-Retention Pipette Tips. 96 Tips/Rack. 10 Racks/Pack. 5 Packs/Case	96	10	5
LTR-20-SL	20uL LTS-Compatible Low-Retention Pipette Tips. Sterile. 96 Tips/Rack. 10 Racks/Pack. 5 Packs/Case	96	10	5
LTR-200-L	200uL LTS-Compatible Low-Retention Pipette Tips. 96 Tips/Rack. 10 Racks/Pack. 5 Packs/Case	96	10	5
LTR-200-SL	200uL LTS-Compatible Low-Retention Pipette Tips. Sterile. 96 Tips/Rack. 10 Racks/Pack. 5 Packs/Case	96	10	5
LTR-1000-L	1000uL LTS-Compatible Low-Retention Pipette Tips. 96 Tips/Rack. 10 Racks/Pack. 5 Packs/Case	96	10	5
LTR-1000-SL	1000uL LTS-Compatible Low-Retention Pipette Tips Sterile. 96 Tips/Rack. 10 Racks/Pack. 5 Packs/Case	96	10	5