

## Oxford Accutip Racked (OAR) Pipette Tips Sterile Low Retention Filter Tips

The Oxford Accupet Racked (OAR) tip line is our newest addition to our universal pipette tip offering. Designed to precisely fit most universal-style single and multichannel pipettes, the OAR line is where form meets function. All Oxford pipette tips are manufactured in a fully automated ISO 9001 facility and certified free of RNase, DNase, PCR-inhibitors, endotoxins, and pyrogens. Sterile tips are certified sterile by gamma irradiation to SAL 10<sup>-6</sup> (ISO11137).

### Features and Benefits

- Ultra-low retention and tapered tips to minimize sample retention and maximize sample delivery
- Universal-friendly to fit virtually all universal-style pipettes
- Ultra-clear and graduated for easy sample visualization
- Manufactured from medical-grade polypropylene in a class 100,000 cleanroom environment under ISO13485 and ISO9001 standards
- Shelf life: 3 years from manufacture date



Cat. No.	Description	Sterile	Qty/Pack	Qty/Box	Qty/Case
OAR-10XL-SLF	10 $\mu$ L, Long, Natural, With filter, Sterilized, Packaged in rack box	Yes	96	960	1920
OAR-20-SLF	20 $\mu$ L, Natural, With filter, Sterilized, Packaged in rack box	Yes	96	960	1920
OAR-100-SLF	100 $\mu$ L, Natural, With filter, Sterilized, Packaged in rack box	Yes	96	960	1920
OAR-200-SLF	200 $\mu$ L, Natural, With filter, Sterilized, Packaged in rack box	Yes	96	960	1920
OAR-300-SLF	300 $\mu$ L, Natural, With filter, Sterilized, Packaged in rack box	Yes	96	960	1920
OAR-1000XL-SLF	1000 $\mu$ L, Long, Natural, With filter, Sterilized, Packaged in rack box	Yes	96	960	1920

