


Maximize the use of precious bench-top space.

**Oxford MHS-10** (High Performance multi-position hotplate stirrer)

## 300 to 1500 rpm with Digital Timer

- Closed design allows for easy cleaning and protects equipment from accidental spills
- Corrosion resistant stainless steel with silicone sheet.
- Timer setting from 1 to 100 minutes and infinite mode
- Anti skid matt for maximum stability
- Digital display for added convenience
- Programmable pulse mode
- 10 stirring positions
- Variable temperature setting from ambient to 120C
- Variable speed settings ranging from 300-1500rpm
- IP 42 class complaint rating for enhanced durability



5 Year Warranty

For more details, visit us at [www.oxfordlp.com](http://www.oxfordlp.com)

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

\*OEM request can be entertained. Specifications can be changed without notice for quality improvement.